



Baby scale "Home"

8310.01

- ▶ Handy and non-tilting baby scale
- ▶ Easy to clean
- ▶ Tare and hold function
- ▶ Energy-saving auto-off function

Features

The Soehnle Professional Baby Scale 8310 convinces with its combination of simple design, easy operation and maintenance-free operation. The scale is ideal for mobile use due to its compactness and low weight.

- ▶ Easy to read LC display, 14 mm digit height
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Easy to clean plastic housing
- ▶ Switchable to kg/lb units
- ▶ Energy-saving auto-off function
- ▶ Battery operation with four AA alkaline batteries (included in scope of delivery)
- ▶ Battery status indicator

Products

Order number Variation	Approvable	Division g	Maximum load kg
8310.01.001 Baby scale "Home"	No	10	20

Accessories

Order number
Variation

Description

5040.02.001



Carry bag black, with carrying strap and durable zipper
