



# Baby scale "Multina Comfort"

8352

- ▶ Detachable baby tray
- ▶ With attachable baby height rod (optional)
- ▶ Simple operation with only one button
- ▶ Scale can be used to weigh children up to 20 kg



## Features

The Soehnle Professional Baby Scale 8352 convinces by its simple but high-class design in combination with highest quality - made in Germany. The intuitive one-button operation, the tare function and the easy-to-clean construction make the baby scale a popular companion for midwives, mothers and birth centres.



The highlight of this scale is the removable baby car seat. This means that even small children up to 20 kg can be weighed on the scale in a standing or sitting position.

- ▶ Simple operation with only one button
- ▶ With detachable baby tray
- ▶ Scale can be used for small children up to 20 kg
- ▶ Easy to read LC-Display, 13 mm digit height
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function: The weight value will be held for a few seconds after unloading
- ▶ Units in kg or lb
- ▶ Energy-saving auto-off function
- ▶ Battery operation with 9 V battery (included in scope of delivery)

## Products

Order number Variation	Approvable	Division g	Maximum load kg
8352.01.001 Baby scale "Multina Comfort"	No	10	20

## Accessories

Order number Variation		Description
5005.11.001		Baby height rod
5040.03.001		Carry bag black, with detachable shoulder strap and durable zipper
5020.01.001		Baby tray (replacement)