

Professional medical devices

FAZZINI



DIATROM®80D

Dim.: 25,4 x 28,8 x 10,4 h cm

EN Electro-surgical unit

Monopolar and bipolar radio-surgical unit suitable for micro-coagulative operations, small and medium surgery. Different functions can be selected from the control panel: CUT, BLEND, COAG FORCED, COAG SOFT, BIPOLAR. Its digital display guarantees reliability and simple use

FR Bistouri électrique

Bistouri électrique monopolaire/bipolaire conçu pour les interventions de micro coagulation, de petite et moyenne chirurgie. Son panneau de commande permet de sélectionner la fonction désirée simplement: CUT, BLEND, COAG FORCED, COAG SOFT, BIPOLAR. L'affichage numérique de la puissance délivrée rend les actes sûrs et reproductibles

IT Elettrobisturi

Elettrobisturi monopolare/bipolare progettato per micro coagulazione, chirurgia piccola e media. Il pannello di controllo consente di selezionare semplicemente la funzione desiderata: CUT, BLEND, COAG FORCED, COAG SOFT, BIPOLAR. Il display digitale della potenza erogata rende gli atti sicuri e riproducibili

Technical specifications / Spécifications techniques / Specifiche tecniche

CODE	10100.10 NEW
Frequency / Fréquence / Frequenza	600 kHz
Max power / Puissance max / Potenza max.	CUT 80 W - 250 Ohm BLEND 60 W - 200 Ohm FORCED COAG 50 W - 150 Ohm SOFT COAG 40 W - 100 Ohm BIPOLAR 30 W - 100 Ohm
Patient circuit / Électrode neutre / Elettrodo neutro	F
Power supply / Alimentation / Alimentazione	115-230Vca / 50-60Hz
Power max absorption / Puissance max. absorbée / Potenza assorbita	200 VA
Weight / Poids / Peso	5 kg

Standard accessories / accessoires / accessori	Qty
00100.03 Cord, 2 m / Câble 2 m / Cavo 2 m	1
F4243 Reusable sterile handpiece / Porte-électrodes monopolaire stérile / Portaelettrodo monopolare sterile	1
00404.08 Cable for cylindrical electrode / Câble de connexion électrode / Cavo di collegamento dell'elettrodo	1
5365 Small neutral plate / Plaque neutre / Piastra neutra	1
00304.00 Pedal switch / Pédale interrupteur / interruttore a pedale impermeabile	1
00500.00 10 Kit of 10 electrodes / Kit de 10 électrodes assorties / Kit di 10 elettrodi assortiti	1
MA280 User's manual / Notice d'emploi / Manuale d'uso	1