



AUDIOMETER 9910



*DD45 headphones +
Peltor® audio cups*

The digital audiometer **9910** allows the healthcare professional to establish an **accurate screening** of the hearing losses.

REVIEW : JANUARY 2021



AUDIOMETER 9910



Device to be put on a desk, works on AC power (or with rechargeable battery as an option), practical and easy to use



On 2 large backlit LCD displays, you can see :
Sound level, frequency, answer of the patient (with the hand switch button)



- Continuous or pulsed sound
- Frequencies from 125 to 8000 Hz, output level range from -10 to 100 HLdB



Delivered into a shock-proof briefcase.
With 2 headphones to choose from :

- A standard headphone or an acoustic highly isolated headphone
- A patient response hand switch



- Device designed and manufactured in France
- Warranty : 3 Years *



End user manual: English, German, Italian, Spanish, Portuguese, Dutch, Hungarian, Romanian, French



Photo with standard headphones

WARRANTY 3 YEARS (EXCEPT FOR THE HEADPHONES) – FRENCH DESIGN AND MANUFACTURING

TECHNICAL INFORMATION

- **Type** : 4 – Screening and Control (according to EN 60645-1) Medical Marking EC 0459 - Class IIa
- **Output** : AC – Air Conduction
- **Frequencies range (11 values)** : 125, 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000, 8000 Hz (+/- 1%)
- **Sound level** : from -10 HLdB to 100 HLdB

- **Increment of sound level** : 5 dB
- **Power supply** : AC power adaptor 230 / 15 V – 500 mA. Integrated rechargeable battery as an option
- **Stimulus** : continuous or pulsed sound
- **Dimensions of the audiometric box** : 210 x 130 x 30/60 mm
- **Packaging** : Shock proof briefcase 280 x 310 x 100 mm
- **Total weight** : 1,5 Kg (500g for the audiometric box only)

Electronica Medical offers a complete range of devices for screening the hearing losses, which are dedicated to the educational and occupational medicine.

Electronica Technologies reserves itself the right to modify the specifications indicated without prior notice. The photos are not contractual.

Headquarters and Industrial Center : Chemin du Pré Long – 03200 Abrest – France. Tél. +33 (0)4 70 96 52 90 - Fax +33 (0)4 70 96 01 23 - contact@electronica-medical.com

* Except for the headphones