

 **semper
med**

derma PF

A highly valued surgical glove

The sempermed® derma PF offers good reasons to switch to powder-free gloves.



Easy to don



Secure grip



Economically-sound solution



Skin-friendly



Easy to don

The glove is easy to put on thanks to the synthetic lining.



Secure grip

Sempermed® derma PF offers a confident grip and thus combines safety and tactile properties.



Economically-sound solution

Allergies and other health risks pose risks to end-users and strain healthcare systems. This powder-free glove eliminates powder-related health issues, reducing long-term costs.





Skin-friendly

This powder-free glove provides improved safety for staff, patients and visitors as there are no powder-related dermatological, allergic and toxic health risks.

Technical data

Type	sterile surgical glove for single use, powder-free, with synthetic lining				
Color	natural white				
Size / Overall length as per EN 455-2	5½, 6 and 6½ / 7, 7½ and 8 / 8½ and 9	median 260 mm / median 270 mm / median 280 mm			
Wall thickness measured in single layer in the palm area	0.16 - 0.21 mm				
Barrier performance	AQL 0.65				
Material	natural rubber latex (NR)				
Glove shape	fully anatomical with rolled rim				
Force at break as per EN 455-2	≥ 12 N, exceeding EN 455-2 (≥ 12 N instead of ≥ 9 N)				
Shelf life in original packaging if stored according to product specification	3 years				
Sterilization	radiation <table border="1"><tr><td>STERILE</td><td>R</td></tr></table>			STERILE	R
STERILE	R				
Labeling & conformity to standards	EN 455-1/-2/-3/-4, EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, ISO 374-5, EN 421, EN 556, ISO 2859-1, ISO 15223-1, ASTM F1671/F1671M				
Purpose	Single-use medical device class IIa according to MD Directive 93/42/EEC				
	Single-use protective glove – PPE category III (protection against chemical substances for limited time)* according to PPE Regulation (EU) 2016/425				





Packaging

- Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in peel pack.
- In dispenser carton with sterilization indicator: 50 pairs
- In transport carton with sterilization indicator: 300 pairs

Article numbers

size 5½	822751522
size 6	822751602
size 6½	822751622
size 7	822751702
size 7½	822751722
size 8	822751802
size 8½	822751822
size 9	822751902